

Die Cast Machinery, LLC Inventory number: DCM-4428

Used Bruker Model: Q6 Columbus QC/V 32 Channel Arc/Spark Optical Emission Spectrometry (OES) analyzer Metal Analytic Spectrometer

Status:	Sold
Manufacturer:	Bruker
Model:	Q6 Columbus QC/V
Size:	32 Channel
Inventory Number:	DCM-4428

Bruker - Q6 Columbus QC/V 32 Channel Arc/Spark Optical Emission Spectrometry (OES) analyzer Metal Analytic Spectrometer, Controls

[View this machine online: DCM-4428](#)



Specifications

	U.S. Standard	Metric
Model:	QC/V	QC/V
Current Baseline Element:	Aluminum (AL)	Aluminum (AL)
Number of Channels:	32	32
Number of Elements Inspected as Set-Up:	13	13
Elements Inspected:	Si, Fe, Cu, Mn, Mg, Ni, Zn, Sn, Cr, Cd, Pb, Ti, Al	Si, Fe, Cu, Mn, Mg, Ni, Zn, Sn, Cr, Cd, Pb, Ti, Al

Features & Accessories:

- Currently Configured for Aluminum Base Materials
- Baseline element: AL (aluminum)
- Number of element inspected: Thirteen (13)
- Elements Inspected: Si, Fe, Cu, Mn, Mg, Ni, Zn, Sn, Cr, Cd, Pb, Ti, Al

Terms of Sale:

Price: [sold](#)

Shipping Terms: Ex-works

Questions? Call (847) 360-9170 and mention DCM-4428.

Photos - DCM-4428

Click any photo to open it full size on our website.




[Photo 1 of DCM-4428 - view online](#)




[Photo 2 of DCM-4428 - view online](#)



[Photo 3 of DCM-4428 - view online](#)



[Photo 4 of DCM-4428 - view online](#)




[Photo 5 of DCM-4428 - view online](#)



[Photo 6 of DCM-4428 - view online](#)

Photos - DCM-4428 (continued)

Click any photo to open it full size on our website.



[Photo 7 of DCM-4428 - view online](#)



[Photo 8 of DCM-4428 - view online](#)




[Photo 9 of DCM-4428 - view online](#)



[Photo 10 of DCM-4428 - view online](#)




[Photo 11 of DCM-4428 - view online](#)



[Photo 12 of DCM-4428 - view online](#)

Photos - DCM-4428 (continued)

Click any photo to open it full size on our website.



[Photo 13 of DCM-4428 - view online](#)



[Photo 14 of DCM-4428 - view online](#)

Description:

Used, Bruker Model Q6 Columbus QC/V 32 Channel Arc/Spark Optical Emission Spectrometry (OES) analyzer Metal Analytic Spectrometer For_Sale, Features Include: Baseline element: AL (aluminum), Number of element inspected: Thirteen (13), Elements Inspected: Si, Fe, Cu, Mn, Mg, Ni, Zn, Sn, Cr, Cd, Pb, Ti, Al, controls, second Hand, DCM-4428

PURCHASER'S (Buyer's) RESPONSIBILITY AND INDEMNITY: Used machines, equipment or any part thereof, may not incorporate approved activating mechanisms, operator safety devices or safety guards as required by OSHA or otherwise. It shall be a Purchaser's (Buyer's) responsibility to ensure that any and all merchandise or equipment purchased from Die Cast Machinery, LLC or its affiliates is installed and operated in a proper and safe manner. Purchaser (Buyer) also acknowledges that it may have to install or change guards, safeties, warnings or other components to ensure that the merchandise or machines purchased will conform to all laws, codes, regulations, ordinances, statutes, insurance requirements and industry standards. Purchaser agrees to defend, indemnify and hold harmless Die Cast Machinery, LLC from and against all suits, claims, costs, damages and expenses, including, but not limited to, reasonable attorneys' fees, arising out of, or in connection with, the transportation, purchase, ownership, or use of the merchandise or equipment sold hereunder. Information has been obtained from sources deemed reliable, buyers should avail themselves of the opportunity to inspect the items before purchase.